


## Three-Prong Cannula Instructions

This cannula uses its own tubing to secure the prongs in place under your nose.




### How to Use Your Cannula

- **Place the cannula under your nose:** Insert the shorter prongs gently into your nostrils, with the longer prong pointing down toward your mouth.
- **Route the tubing:** Bring the two individual tubes up and over the tops of your ears, just like you would wear eyeglasses.
- **Secure under your chin:** Guide the tubing ends down to meet under your chin.
- **Adjust for comfort:** Adjust the slider/clip to tighten or loosen the tubing until the cannula feels snug and comfortable.

## Confirm Airflow is Detected

Before you begin the recording, speak or blow into the nasal cannula. The LED light on the airflow sensor (microphone) will turn green, indicating proper placement and function.



### Tips for Your Comfort

- The nasal prongs should rest comfortably in your nostrils without causing pressure. You can trim the tips if they feel too long.
- For further stability, medical tape may be used to secure the cannula to your cheeks if needed.

If you have any questions, please email **Support@SnapDiagnostics.com** or call us at **(847) 777-0000**.